



SURGICAL SUTURES



THREAD TYPE

Polyglycolic acid (PGA)
 Chromic catgut (CC)
 Plain catgut (CP)
 Polypropylene (PM)
 Polyester braided (PB)
 Silk braided (SK)
 Nylon (NL)
 Stainless steel wire (SW)
 Ligature (LR)

THREAD LENGTH

45cm, 60cm, 75cm, 100cm, 125cm, 150cm

THREAD DIAMETER

8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2, 3

NEEDLE LENGTH







6mm, 8mm, 12mm, 18mm, 22mm,
 30mm, 35mm, 40mm, 50mm

NEEDLE CURVATURE


Straight, 1/2 circle, 1/2 circle (double),
 1/4 circle, 1/4 circle (double)
 3/8 circle, 3/8 circle (double), 5/8 circle,
 loop round

CROSS SECTION

Round bodied, round bodied (heavy),
 curved cutting, curved cutting (heavy),
 Reverse cutting, reverse cutting (heavy),
 tapercut, micro-point spatula curved

POINT TYPE	SYMBOL
 TAPER POINT	●
 BLUNT TAPER POINT	○
 CUTTING EDGE	▲
 REVERSE CUTTING EDGE	▼
 TAPERCUT	⊕
 MICRO-POINT SPATULA CURVED	◄

SHAPES



SUTURE MATERIALS

