

Laboratory Equipment



Medline Vortex Mixer


- Strong and silent motor
- Touch activation or continuous modes
- Suction cups on base minimise vibration and movement
- Supplied with tube adapter and flask accessories
- Optional accessories are also available to expand the usage of the mixer

Product reference VM-3000MD Specifications

Speed	20 to 3000rpm
Dimensions	130(W)x155(D)x155(H)mm
Weight	3.1 Kg
Power:	220V 50/60Hz



Microcentrifuge tube accessory
VM-3000MT



Microtitre plate accessory
VM-3000MP



Flask accessory
VM-3000FP



Ultrasonic Cleaners

An industrial BLT trembler adapted to produce highly efficient pulsation waves, this range enables powerful and reliable cleaning and the output can be set to high, medium or low in 3 easy stages. The external casing and cover are fully waterproof and corrosion-resistant.

Model	UC-02	UC-10	UC-20
Product reference	BDH53015	BDH53045	BDH53065
Bath volume litres	2.0	10.0	20.0
Temperature range	Amb. ~ to 50°C	Amb. ~ to 70°C	Amb. ~ to 70°C
Ultrasonic power (W)	70	200	400
Dimensions			
Internal (W x L x H)mm	142 x 132 x 104	300 x 240 x 150	495 x 292 x 150
Overall (W x L x H)mm	290 x 210 x 200	49 x 368 x 297	721 x 451 x 297
Net weight kg	2.0	11.0	17.0
Electrical Requirements (230V, 50/60Hz)	0.3A	0.9A	1.7A

Multiple frequency on 3 stages is available for both gentle and intensive cleaning on the UC-10 and UC-20 models.

The Ultrasonic power (W) is increased to the following:

UC-10P - 300

UC-20P - 500



Aspirators

The VE-11 electric aspirator is designed especially for low-vacuum applications and is perfect for filtration.

- Durable & tough design
- Quiet operation level
- Excellent chemical resistance
- Re-circulating water function

Product reference	71015
Pump Capacity	18 L / min
Bath Capacity	9.5 L
Dimensions	
Internal (W x D x H)mm	298 x 210 x 227
Overall (W x D x H)mm	330 x 265 x 390