


**REMI PMDC LAB STIRRER  
RQ-121/D**

**R 2250 /EACH**



**PMDC LAB STIRRERS**

The new PMDC Stirrers are specially designed for a broad spectrum of applications. These models are designed for routine laboratory work, research and development, QA analysis and sample preparations. Stirrers have uniform speed control with digital RPM indicator which gives speed read outs, valuable for applications where process validation and reproducibility is required.

RQ-121/D, RQ-126/D and RQ-129/D models are ideal for day to day laboratory use with aqueous & low viscosity liquids such as light oils.

RQG-121/D, RQG-126/D & RQG-129/D are geared models with hollow shafts in output spindle, permitting length adjustment with ability to mix high viscosity liquids such as heavy oils, lubricants & paints.

**REMI REVOLUTIONARY MICRO CENTRIFUGE  
RM-12 C DX**

**R 4999.00 EACH  
PLUS ROTOR**

**REVOLUTIONARY MICRO  
CENTRIFUGES**

These Centrifuges are ideal for growing RIA, routine application in biochemical and clinical labs, for Hematocrite, Corpuscle Percentage contents in blood, serum analysis and precipitate separation etc. RM-12 is the basic unit with speed regulator. RM-12C has digital speed indicator & 0-15 minutes digital countdown timer.

RM-12C DX has additional features of dynamic brake for quick deceleration & imbalance detector with cutoff in case of uneven loads. RM-12C BL has additional features of brushless induction motor with frequency inverter making pre-setting of speeds possible RM-12C DX & RM-12C BL meets the requirements of IEC-61010 and is available with CE mark, (Conformity to European safety standards)




**B & M SCIENTIFIC cc**  
**Suppliers of Quality Laboratory**

**Equipment & Chemicals**

**Tel : +27 21 933 8132**

**Fax : +27 21 933 8931**

**sales@bmscientific.co.za**



**UNIRACK  
POLYPROPYLENE  
DOUBLE SIDED RACK**

**R 49.00 EACH**



**MEDIA STORAGE BOTTLE  
NON AUTOCLAVABLE  
1000 ml**

**R34.00 EACH**

**PCR RACK  
POLYPROPYLENE**

**R 45.00 EACH**



**YELLOW TIPS  
VOL. 200ul  
GILSON/UNIVERSAL TYPE**

**R 39.00/1000**

**BLUE TIPS  
VOL.1000ul  
GILSON/UNIVERSAL TYPE**

**R 59.00/1000**

**CRYO BOXES  
CARDBOARD  
81 PLACE**

**R 26.00/ EACH**

**CRYO BOXES  
POLYCARBONATE  
81 PLACE  
1ml OR 1.8 CRYO VIALS**

**R 59.00/EACH**

**CENTRIFUGE TUBE**

**VOL. 50 ml, POLYPROPLENE,  
SKIRTED/ CONICAL  
PRINTED GRADUATION,STERILE**

**R 899.00 / 500**

**MICROSCOPE COVER GLASS**

**18 X 18 mm R 8.00/100  
22 X 40 mm R 12.00 /100**

**MICROSCOPE SLIDE FROSTED**

**R 8.50 /50**

**LENS PAPER  
10 X 15 cm**

**R 14.00/ BOOKLET 100 LVS**

**TEST TUBE  
VOL. 5 ml  
12 X 75mm**

**R 99.00 / 1000 POLYSTYRENE  
R 120 / 250 GLASS**

**SCHOTT EQUIVALENT MEDIA BOTTLE  
AUTOCLAVABLE**

**1000ml R 39.90 /EACH  
500ml R 28.00 /EACH  
250ml R 21.00 /EACH**

**PASTEUR PIPETTE  
GLASS  
23 cm**

**R 89.90 / 250**

**MICROCENTRIFUGE TUBES**

**AUTOCLAVABLE, NON STERILE  
HINGED LID & WRITING PANEL-  
GRADUATED**

**1.5 ml R 89.00 / 1000  
0.6 ml R 90.00 / 1000  
0.2 ml R 79.00 / 1000**

**NON STERILE SYRINGE FILTERS  
0.22 um**

**R 89.90 / 50**

**SEROLOGICAL PIPETTES  
5 ml**

**R 929.00 / 500**

**SWAB COTTON TIP IN TUBE  
150mm  
STERILE**

**R 149 / 100**

**BEAKERS  
BOROSILICATE GLASS  
LOW FORM  
100 ml**

**R 6.00 / EACH**

**WE ARE SUPPLIERS OF LABORATORY AND RESEARCH CHEMICALS,  
EQUIPMENT, CONSUMABLES AND MEDICAL SUPPLIES**

**FOR MORE INFORMATION ABOUT OUR SPECIALS, PRODUCTS OR GENERAL QUES-  
TIONS, PLEASE DO NOT HESITATE TO CONTACT US, WE WILL GLADLY ASSIST**

**ALL ITEMS ARE WHILE STOCKS LAST  
PLEASE CONFIRM STOCK AVAILIBLTY**

**REGARDS THE B & M SCIENTIFIC TEAM**